prosthetic.line

CM Connect

COBALT-CHROME BASE.

Instructions for use.

CHROME-COBALT BASIS.

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BASE COBALT-CHROME.

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BASE COBALTO-CROMO.

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
BASE DE CROMO-COBALTO.

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CoCr Base

Extended instructions for use.

A prefabricated CoCr base with plastic burnout sleeve is the best alternative for precious alloys. Thanks to serial machine milling the implant interface creates a perfect surface which is important for successful implant restoration. It is crucial for implant screw durability and for any bacterial leakage along the implant-abutment connection. The photo in Scanning Electron Microscopy (SEM) shows the difference between adaptation of two connections: prefabricated (left) and after laboratory casting (right). [C. Domingues da Silveira].

Depending on implant brand and connection type you can choose from a range of our products:


- can be used in bone level restorations or from abutment level (Multiunit),
- Engaging connection for custom abutments, telescopic crowns or single restorations fixed with occlusal screw,
- Non-engaging version for bridges fixed with occlusal screws or mesostructure in combined restorations.

They can be used for over-casting technique as well as for laser welding. The castable plastic sleeve is made of sanitary polyoxymethylene, which burns out without any residues.

MODELLING TECHNIQUE

1. After choosing the correct interface the piece is placed on the laboratory analog and fixed with prosthetic screw by hand.
2. The plastic burn out sleeve is trimmed according to the anatomical and occlusal requirements.
3. Wax-up for ceramic is made following the conventional crown and bridge requirements.
4. Create a collar of at least 0.3 mm thickness above the CoCr base (as seen on image). It prevents cracks in ceramic layer.
5. The wax-up is removed and before placing casting sprues the CoCr base should be cleaned with a cotton swab and alcohol.
6. Casting sprues should be placed possibly in the long axis of the crown to avoid the bubble concentration inside the construction. The investment should freely pass through the screw channel. The CoCr base needs to be placed away from casting ring heating center. The use of speed investment is not recommended.
7. The investment expansion should be adjusted to zero expansion. Heat parameters need to be tuned as for conventional CoCr casting to avoid any imperfections in the construction.
8. It is better to deflask the casting mould using an acid bath in an ultrasonic cleaner. Sandblasting can damage the settlement.
9. During neck polishing the connection surface needs to be covered with analog. It prevents any damage of the settlement part.


CM CONNECT CoCr INTERFACE PHYSICAL AND CHEMICAL PROPERTIES

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