

CM Connect

Library guide for hybrid abutments,
bridges, bars and prosthesis.

exocad

Exocad library for implant level scan and abutment level scan























For this CAD system the library draws distinction between the following groups:


exocad


- CM Connect Ti-base – for restorations connected with the implant via the Ti-base . The scan needs to be performed with the Scan Body for the specific implant connection. **Note:** This group includes parts **scanned from abutment level** - Multi-Unit abutment, UniAbutment® and synOcta® abutment - with its specific Scan Abutment (marked almond in the table).
- CM Connect MULTI (or MU) - for restorations connected directly with Multi-Unit transepithelial abutment. The scanning needs to be performed with the Scan Body for the specific implant connection. The desired Multi-Unit abutment will be mounted virtually in the CAD System. Please note that in this case the closing line between the prosthesis and the transepithelial abutment depends on the proper kit selection (abutment height in mm).
- CM Connect MU-P (or MULTI PLUS) Ti-base - for restorations connected with the Multi-Unit transepithelial abutment via the Ti-base . The scanning needs to be performed with the Scan Body for the specific implant connection. The desired Multi-Unit abutment will be mounted virtually in the CAD System. Please note that in this case the closing line between the prosthesis and the transepithelial abutment depends on the proper kit selection (abutment height in mm).
- CM UNIABT-P Ti-base - for restorations connected with the UniAbutment® via Ti-base. The scanning needs to be performed with the Scan Body for ASTRA TECH™ implant connection. The desired UniAbutment® will be mounted virtually in the CAD System. Please note that in this case the closing line between the prosthesis and the transepithelial abutment depends on the proper kit selection (abutment height in mm).
- CM AS UNIABT - for restorations connected directly with the UniAbutment®. The scanning needs to be performed with the Scan Body for ASTRA TECH™ implant connection. The desired UniAbutment® will be mounted virtually in the CAD System. Please note that in this case the closing line between the prosthesis and the transepithelial abutment depends on the proper kit selection (abutment height in mm).
- CM Connect AURUMBase® – for restorations connected with the implant via the AURUMBase® (see Table 1.4). The scanning needs to be performed with the Scan Body for the specific implant connection. **Note:** This group includes parts **scanned from abutment level** - Multi-Unit abutment - with its specific Scan Abutment (marked almond colour in the table).
- CM Connect Pre-milled Ti – for restorations connected with the implant via the pre-milled abutment made in Titanium (see Table 1.6). The scanning needs to be performed with the Scan Body for the specific implant connection. **Note:** The maximum width of the piece is restricted to 10 mm.
- CM Connect Pre-milled CoCr – for restorations connected with the implant via the pre-milled abutment made in CoCr (see Table 1.7). The scanning needs to be performed with the Scan Body for the specific implant connection. **Note:** The maximum width of the piece is restricted to 14 mm.

Take care to select the correct CM Scan Abutment and the corresponding library file, as otherwise the milled prosthesis can differ from desired outcome and errors can occur, for example the level of interface settlement.


Table 1.1
CM Connect TI-Base


Name Library	Subtype	Implant Type	Analogue	Scan Abutments	
Nobel Active®		3.0	41.14.068	41.52.068	
				41.53.100	
					41.52.068
					41.52.200
		NP/3.5	41.14.041		41.50.041 41.52.041
					41.53.100
					41.50.041 41.52.041
					41.52.200
		RP/4.3	41.14.042		41.50.042 41.52.042
					41.53.100
					41.50.042 41.52.042
					41.52.200
	WP/5.5	41.14.069		41.52.069	
				41.53.100	
				41.52.069	
				41.52.200	
Nobel Replace Select™		NP/3.5	41.14.004	41.50.004 41.52.004	
				41.53.100	
					41.50.004 41.52.004
					41.52.200
		RP/4.3	41.14.005		41.50.005 41.52.005
					41.53.100
					41.50.005 41.52.005
					41.52.200
		WP/5.0	41.14.006		41.50.006 41.52.006
					41.53.100
					41.50.006 41.52.006
					41.52.200
	6.0	41.14.097		41.52.097	
				41.53.100	
				41.52.097	
				41.52.200	
Nobel Brånemark®		NP/3.5	41.14.001	41.50.001 41.52.001	
				41.53.100	
					41.50.001 41.52.001
					41.52.200
		RP/4.1	41.14.002		41.50.002 41.52.002
					41.53.100
					41.50.002 41.52.002
					41.52.200
		WP/5.0	41.14.003		41.50.003
					41.52.003
					41.50.003
					41.52.003
























 non-engaging


 non-engaging


 Scan Abutment for desktop scanning
41.50.0XX

 Scan Abutment for intra-oral scanning
41.52.0XX


 Scan Abutment on Ti base
41.53.100 41.52.200


Name Library	Subtype	Implant Type	Analogue	Scan Abutments	
3i Osseotite®		NP/3.4	41.14.011	41.50.011 41.52.011 41.53.100	
					41.50.011 41.52.011 41.52.200
		RP/4.1	41.14.012	41.50.012 41.52.012 41.53.100	
					41.50.012 41.52.012 41.52.200
		WP/5.0	41.14.013	41.50.013 41.52.013 41.53.100	
					41.50.013 41.52.013 41.52.200
3i Certain®		NP/3.4	41.14.014	41.50.014 41.52.014 41.53.100	
					41.50.014 41.52.014 41.52.200
		RP/4.1	41.14.015	41.50.015 41.52.015 41.53.100	
					41.50.015 41.52.015 41.52.200
		WP/5.0	41.14.016	41.50.016 41.52.016 41.53.100	
					41.50.016 41.52.016 41.52.200
Straumann® Soft Tissue Level		NNC/3.5	41.14.070	41.52.069 41.53.100	
					41.52.069 41.52.200
					
		RN/4.8	41.14.009	41.50.009 41.52.009	
					41.50.010 41.52.010
					
		WN/6.5	41.14.045	41.50.045 41.52.045	
					41.50.045 41.52.045
					
		WN on synOCTA®	41.14.046	41.50.015 41.52.015	
					41.50.015 41.52.015

 non-engaging


 non-engaging

 Scan Abutment
for desktop scanning
41.50.0XX


 Scan Abutment
for intra-oral scanning
41.52.0XX

 Scan Abutment
on Ti base
41.53.100 41.52.200

Name Library	Subtype	Implant Type	Analogue	Scan Abutments
Straumann® Bone Level	CM - STR BL Ti-Base (Desktop or Intraoral or on Ti base)	NC/3.3	41.14.043	41.50.043 41.52.043 41.53.100
				41.50.043 41.52.043 41.52.200
		RC/4.1	41.14.044	41.50.044 41.52.044 41.53.100
				41.50.044 41.52.044 41.52.200
		SRA	41.14.098	41.52.098
				41.52.098
Zimmer Screw-Vent®	CM - ZI SCV Ti-Base (Desktop or Intraoral or on Ti base)	NP/3.5	41.14.017	41.50.017 41.52.017 41.53.100
				41.50.017 41.52.017 41.52.200
		RP/4.5	41.14.018	41.50.018 41.52.018 41.53.100
				41.50.018 41.52.018 41.52.200
		WP/5.7	41.14.019	41.50.019 41.52.019 41.53.100
				41.50.019 41.52.019 41.52.200
Astra tech implant system™ EV	CM - AS EV Ti-Base (Desktop or Intraoral or on Ti base)	EV/3.6	41.14.059	41.50.059 41.52.059 41.53.100
				41.50.059 41.52.059 41.52.200
		EV/4.2	41.14.060	41.50.060 41.52.060 41.53.100
				41.50.060 41.52.060 41.52.200
		EV/4.8	41.14.061	41.50.061 41.52.061 41.53.100
				41.50.061 41.52.061 41.52.200
Uniabutment®	41.14.058	41.50.058		
		41.50.058 41.52.058		




non-engaging



non-engaging






















Scan Abutment
for desktop scanning
41.50.OXX





Scan Abutment
for intra-oral scanning
41.52.OXX




Scan Abutment
on Ti base
41.53.100 41.52.200


Name Library	Subtype	Implant Type	Analogue	Scan Abutments	
Astra tech osseospeed™		NP/3.0	41.14.023	41.52.023	
				41.53.100	
					41.50.059 41.52.059
					41.52.200
		RP/3.5-4.0	41.14.024		41.50.024 41.52.024
					41.53.100
					41.50.024 41.52.024
					41.52.200
		WP/4.5-5.0	41.14.025		41.50.025 41.52.025
					41.53.100
				41.50.025 41.52.025	
				41.52.200	
	Uniabutment® 20°		41.14.036	41.50.036	
				41.52.036	
	Uniabutment® 45°		41.14.037	41.50.037	
				41.52.037	
Friadent XIVE®		NP/3.4	41.14.038	41.50.038 41.52.038	
				41.53.100	
					41.50.038 41.52.038
					41.52.200
		RP/3.8	41.14.039		41.50.039 41.52.039
					41.53.100
					41.50.039 41.52.039
					41.52.200
		WP/4.5	41.14.040		41.50.040 41.52.040
					41.53.100
				41.50.040 41.52.040	
				41.52.200	
	5.5	41.14.071		41.52.071	
				41.53.100	
				41.52.071	
				41.52.200	
Ankylos® C/X		All platforms, height GH 1,5	41.14.048 41.14.049	41.50.049 41.52.049	
				41.53.100	
					41.50.049 41.52.049
					41.52.200
		All platforms, height GH 2,5	41.14.048 41.14.049		41.50.049 41.52.049
41.53.100					
				41.50.049 41.52.049	
				41.52.200	













 non-engaging


 non-engaging


 Scan Abutment
for desktop scanning
41.50.0XX

 Scan Abutment
for intra-oral scanning
41.52.0XX


 Scan Abutment
on Ti base
41.53.100 41.52.200


Name Library	Subtype	Implant Type	Analogue	Scan Abutments	
Osstem® TS		Mini Yellow	41.14.072	41.50.072 41.52.072	
				41.53.100	
					41.50.072 41.52.072
					41.52.200
		STD Regular	41.14.073		41.50.073 41.52.073
					41.53.100
				41.50.073 41.52.073	
				41.52.200	
Camlog®		3.3	41.14.064	41.52.064	
				41.53.100	
					41.52.064
					41.52.200
		3.8	41.14.065		41.50.065 41.52.065
					41.53.100
					41.50.065 41.52.065
					41.52.200
		4.3	41.14.066		41.50.066 41.52.066
					41.53.100
				41.50.066 41.52.066	
				41.52.200	
	5.0	41.14.067		41.52.067	
				41.53.100	
				41.52.067	
				41.52.200	

















 non-engaging


 non-engaging


 Scan Abutment for desktop scanning
41.50.0XX

 Scan Abutment for intra-oral scanning
41.52.0XX


 Scan Abutment on Ti base
41.53.100 41.52.200

Name Library	Subtype	Implant Type	Analogue	Scan Abutments	
Biohorizns® External		NP/3.7	41.14.020	41.50.020 41.52.020	
				41.53.100	
					41.50.020 41.52.020
					41.52.200
		RP/4.2	41.14.021		41.50.021 41.52.021
					41.53.100
				41.50.021 41.52.021	
				41.52.200	
	WP/5.2	41.14.022		41.50.022	
				41.52.022	
				41.50.022	
				41.52.022	
MIS® Seven		NP/3.7	41.14.056	41.50.056 41.52.056	
				41.53.100	
					41.50.056 41.52.056
					41.52.200
		SP/4.2	41.14.095		41.50.017 41.52.017
					41.53.100
				41.50.017 41.52.017	
				41.52.200	
	WP/5.2	41.14.096		41.50.018 41.52.018	
				41.53.100	
				41.50.018 41.52.018	
				41.52.200	
Megagen AntRidge®			41.14.057	41.50.057	
				41.53.100	
				41.50.057	
				41.52.200	
Multi-Unit		RP	41.14.007	41.50.007	
		41.15.007		41.52.007	
		WP	41.14.047	41.50.047	
		41.15.054		41.52.047	

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
41.50.0XX

 Scan Abutment
for intra-oral scanning
41.52.0XX



 Scan Abutment
on Ti base
41.53.100 41.52.200

Table 1.2
CM Connect Multi-Unit &
Uniabutment plus Ti-Base

	Name Library	Subtype	Platform	Analogue	Scan Abutments
Multi-Unit on Nobel Active™	CM - NB AC MULTI (Desktop or Intraoral)	NP 41.15.007+41.40.007/1 (1,5mm)	NP	41.14.041	41.50.041
		NP 41.15.007+41.40.007/2 (2,5mm)			41.52.041
		NP 41.15.007+41.40.007/3 (3,5mm)			
		RP 41.15.007+41.40.008/1 (1,5mm)	RP	41.14.042	41.50.042
		RP 41.15.007+41.40.008/2 (2,5mm)			41.52.042
RP 41.15.007+41.40.008/3 (3,5mm)					
RP 41.15.007+41.40.008/4 (4,5mm)					
Multi-Unit on Nobel Replace®	CM - NB RP MULTI (Desktop or Intraoral)	NP 41.15.007+41.40.004/1 (1mm)	NP	41.14.004	41.50.004
		NP 41.15.007+41.40.004/2 (2mm)			41.52.004
		NP 41.15.007+41.40.004/3 (3mm)			
		RP 41.15.007+41.40.005/2 (2mm)	RP	41.14.005	41.50.005
		RP 41.15.007+41.40.005/3 (3mm)			41.52.005
	RP 41.15.007+41.40.005/4 (4mm)				
	RP 41.15.007+41.40.005/5 (5mm)				
WP 41.15.054+41.40.006/1 (1mm)	WP	41.14.063	41.50.006		
WP 41.15.054+41.40.006/2 (2mm)			41.52.006		
WP 41.15.054+41.40.006/3 (3mm)					
Multi-Unit on Brånemark System®	CM - NB BM MULTI (Desktop or Intraoral)	NP 41.15.007+41.40.001/1 (1mm)	NP	41.14.001	41.50.001
		NP 41.15.007+41.40.001/2 (2mm)			41.52.001
		NP 41.15.007+41.40.001/3 (3mm)			
		RP 41.15.007+41.40.002/2 (2mm)	RP	41.14.002	41.50.002
		RP 41.15.007+41.40.002/3 (3mm)			41.52.002
	RP 41.15.007+41.40.002/4 (4mm)				
	RP 41.15.007+41.40.002/5 (5mm)				
WP 41.15.054+41.40.003/1 (1mm)	WP	41.14.003	41.50.003		
WP 41.15.054+41.40.003/2 (2mm)			41.52.003		
WP 41.15.054+41.40.003/3 (3mm)					
Multi-Unit on 3i Osseotite®	CM - 3i OS MULTI (Desktop or Intraoral)	NP 41.15.007+41.40.009/1 (1mm)	NP	41.14.011	41.50.011
		NP 41.15.007+41.40.009/2 (2mm)			41.52.011
		NP 41.15.007+41.40.009/3 (3mm)			
Multi-Unit on Bone Level™	CM - STR BL MULTI (Desktop or Intraoral)	NP 41.15.007+41.40.043/1 (1,5mm)	NP	41.14.043	41.50.043
		NP 41.15.007+41.40.043/2 (2,5mm)			41.52.043
		NP 41.15.007+41.40.043/3 (3,5mm)			
		RP 41.15.007+41.40.044/1 (1,5mm)	RP	41.14.044	41.50.044
		RP 41.15.007+41.40.044/2 (2,5mm)			41.52.044
RP 41.15.007+41.40.044/3 (3,5mm)					
RP 41.15.007+41.40.044/4 (4,5mm)					
Multi-Unit on Zimmer Screw-Vent®	CM - ZI SCV MULTI (Desktop or Intraoral)	NP 41.15.007+41.40.019/1 (1mm)	NP	41.14.017	41.50.017
		NP 41.15.007+41.40.019/2 (2mm)			41.52.017
		NP 41.15.007+41.40.019/3 (3mm)			
		NP 41.15.007+41.40.019/4 (4mm)	RP	41.14.018	41.50.018
		RP 41.15.007+41.40.020/1 (1mm)			41.52.018
RP 41.15.007+41.40.020/2 (2mm)					
		RP 41.15.007+41.40.020/3 (3mm)			



Name Library	Subtype	Platform	Analogue	Scan Abutments		
Uniabutment® on Astra Tech™ Implant System™ EV	CM - AS EV MULTI (Desktop or Intraoral)	EV3.6 41.15.058+41.40.059/1 (1mm)	41.14.059	41.50.059		
		EV3.6 41.15.058+41.40.059/2 (2mm)		41.52.059		
		EV3.6 41.15.058+41.40.059/3 (3mm)				
		EV3.6 41.15.058+41.40.059/5 (5mm)				
		EV4.2 41.15.058+41.40.060/1 (1mm)	EV 4.2 Uniaabutment®	41.14.060	41.50.060	
					EV4.2 41.15.058+41.40.060/2 (2mm)	41.52.060
					EV4.2 41.15.058+41.40.060/3 (3mm)	
					EV4.8 41.15.058+41.40.061/1 (1mm)	EV 4.8 Uniaabutment®
		EV4.8 41.15.058+41.40.061/2 (2mm)	41.52.061			
		EV4.8 41.15.058+41.40.061/3 (3mm)				
EV4.8 41.15.058+41.40.061/5 (5mm)						
Uniabutment® on Astra Tech™ Osseospeed®	CM - AS MULTI (Desktop or Intraoral)	RP 20 41.15.036+41.40.013/0 (0,5mm)	41.14.024	41.50.024		
		RP 20 41.15.036+41.40.013/1 (1mm)		41.52.024		
		RP 20 41.15.036+41.40.013/2 (2mm)				
		RP 20 41.15.036+41.40.013/4 (4mm)				
		RP 20 41.15.036+41.40.013/6 (6mm)				
		WP 20 41.15.036+41.40.015/0 (0,5mm)	WP Uniaabutment® 20°	41.14.025	41.50.025	
		WP 20 41.15.036+41.40.015/1 (1mm)			41.52.025	
		WP 20 41.15.036+41.40.015/2 (2mm)				
		WP 20 41.15.036+41.40.015/4 (4mm)				
		WP 20 41.15.036+41.40.015/6 (6mm)				
		RP 45 41.15.037+41.40.014/0 (0,5mm)	RP Uniaabutment® 45°	41.14.024	41.50.024	
		RP 45 41.15.037+41.40.014/1 (1mm)			41.52.024	
		RP 45 41.15.037+41.40.014/2 (2mm)				
		RP 45 41.15.037+41.40.014/4 (4mm)				
		RP 45 41.15.037+41.40.014/6 (6mm)				
		WP 45 41.15.037+41.40.016/0 (0,5mm)	WP Uniaabutment® 45°	41.14.025	41.50.025	
WP 45 41.15.037+41.40.016/1 (1mm)	41.52.025					
WP 45 41.15.037+41.40.016/2 (2mm)						
WP 45 41.15.037+41.40.016/4 (4mm)						
WP 45 41.15.037+41.40.016/6 (6mm)						






Table 1.3
CM Connect Multi-Unit &
Uniabutment direct

	Name Library	Subtype	Platform	Analogue	Scan Abutments
Multi-Unit on Nobel Active™	CM - NB AC MULTI (Desktop or Intraoral)	NP 41.40.007/1 (1,5mm) NP 41.40.007/2 (2,5mm) NP 41.40.007/3 (3,5mm)	NP	41.14.041	41.50.041 41.52.041
		RP 41.40.008/1 (1,5mm) RP 41.40.008/2 (2,5mm) RP 41.40.008/3 (3,5mm) RP 41.40.008/4 (4,5mm)	RP	41.14.042	41.50.042 41.52.042
Multi-Unit on Nobel Replace®	CM - NB RP MULTI (Desktop or Intraoral)	NP 41.40.004/1 (1mm) NP 41.40.004/2 (2mm) NP 41.40.004/3 (3mm)	NP	41.14.004	41.50.004 41.52.004
		RP 41.40.005/2 (2mm) RP 41.40.005/3 (3mm) RP 41.40.005/4 (4mm) RP 41.40.005/5 (5mm)	RP	41.14.005	41.50.005 41.52.005
		WP 41.40.006/1 (1mm) WP 41.40.006/2 (2mm) WP 41.40.006/3 (3mm)	WP	41.14.0063	41.50.006 41.52.006
Multi-Unit on Brånemark System®	CM - NB BM MULTI (Desktop or Intraoral)	NP 41.40.001/1 (1mm) NP 41.40.001/2 (2mm) NP 41.40.001/3 (3mm)	NP	41.14.001	41.50.001 41.52.001
		RP 41.40.002/2 (2mm) RP 41.40.002/3 (3mm) RP 41.40.002/4 (4mm) RP 41.40.002/5 (5mm)	RP	41.14.002	41.50.002 41.52.002
		WP 41.40.003/1 (1mm) WP 41.40.003/2 (2mm) WP 41.40.003/3 (3mm)	WP	41.14.003	41.50.003 41.52.003
Multi-Unit on 3i Osseotite®	CM - 3i OS MULTI (Desktop or Intraoral)	NP 41.40.009/1 (1mm) NP 41.40.009/2 (2mm) NP 41.40.009/3 (3mm)	NP	41.14.011	41.50.011 41.52.011
Multi-Unit on Bone Level™	CM - STR BL MULTI (Desktop or Intraoral)	NP 41.40.043/1 (1,5mm) NP 41.40.043/2 (2,5mm) NP 41.40.043/3 (3,5mm)	NP	41.14.043	41.50.043 41.52.043
		RP 41.40.044/1 (1,5mm) RP 41.40.044/2 (2,5mm) RP 41.40.044/3 (3,5mm) RP 41.40.044/4 (4,5mm)	RP	41.14.044	41.50.044 41.52.044
Multi-Unit on Zimmer Screw-Vent®	CM - ZI SCV MULTI (Desktop or Intraoral)	NP 41.40.019/1 (1mm) NP 41.40.019/2 (2mm) NP 41.40.019/3 (3mm) NP 41.40.019/4 (4mm)	NP	41.14.017	41.50.017 41.52.017
		RP 41.40.020/1 (1mm) RP 41.40.020/2 (2mm) RP 41.40.020/3 (3mm)	RP	41.14.018	41.50.018 41.52.018




non-engaging



non-engaging



Scan Abutment
for desktop scanning
41.50.0XX



Scan Abutment
for intra-oral scanning
41.52.0XX

	Name Library	Subtype	Platform	Analogue	Scan Abutments
Uniabutment® on Astra Tech™ Implant System™ EV	CM - AS EV MULTI (Desktop or Intraoral)	EV3.6 41.40.059/1 (1mm)	RP Uniabutment®	41.14.059	41.50.059
		EV3.6 41.40.059/2 (2mm)			41.52.059
		EV3.6 41.40.059/3 (3mm)			
		EV3.6 41.40.059/5 (5mm)			
		EV4.2 41.40.060/1 (1mm)	EV 4.2 Uniabutment®	41.14.060	41.50.060
		EV4.2 41.40.060/2 (2mm)			41.52.060
		EV4.2 41.40.060/3 (3mm)			
		EV4.8 41.40.061/1 (1mm)	EV 4.8 Uniabutment®	41.14.061	41.50.061
		EV4.8 41.40.061/2 (2mm)			41.52.061
		EV4.8 41.40.061/3 (3mm)			
		EV4.8 41.40.061/5 (5mm)			
		Uniabutment® on Astra Tech™ Osseospeed®	CM - AS MULTI (Desktop or Intraoral)	RP 20 41.40.013/0 (0,5mm)	RP Uniabutment® 20°
RP 20 41.40.013/1 (1mm)	41.52.024				
RP 20 41.40.013/2 (2mm)					
RP 20 41.40.013/4 (4mm)					
RP 20 41.40.013/6 (6mm)					
WP 20 41.40.015/0 (0,5mm)	WP Uniabutment® 20°			41.14.025	41.50.025
WP 20 41.40.015/1 (1mm)					41.52.025
WP 20 41.40.015/2 (2mm)					
WP 20 41.40.015/4 (4mm)					
WP 20 41.40.015/6 (6mm)					
RP 45 41.40.014/0 (0,5mm)	RP Uniabutment® 45°			41.14.024	41.50.024
RP 45 41.40.014/1 (1mm)					41.52.024
RP 45 41.40.014/2 (2mm)					
RP 45 41.40.014/4 (4mm)					
RP 45 41.40.014/6 (6mm)					
WP 45 41.40.016/0 (0,5mm)	WP Uniabutment® 45°			41.14.025	41.50.025
WP 45 41.40.016/1 (1mm)					41.52.025
WP 45 41.40.016/2 (2mm)					
WP 45 41.40.016/4 (4mm)					
WP 45 41.40.016/6 (6mm)					





























Table 1.4
CM Connect AURUMBase®





















	Name Library	Subtype	Platform	Analogue	Scan Abutments
Nobel Active®	CM - NB AC AURUM BASE (Desktop or Intraoral or on Ti base)		NP	41.14.041	41.50.041 41.52.041 41.53.100
					41.50.041 41.52.041 41.52.200
			RP	41.14.042	41.50.042 41.52.042 41.53.100
					41.50.042 41.52.042 41.52.200
Nobel Replace Select™	CM - NB RP AURUM BASE (Desktop or Intraoral or on Ti base)		NP	41.14.004	41.50.004 41.52.004 41.53.100
					41.50.004 41.52.004 41.52.200
			RP	41.14.005	41.50.005 41.52.005 41.53.100
					41.50.005 41.52.005 41.52.200
			WP	41.14.006	41.50.006 41.52.006 41.53.100
					41.50.006 41.52.006 41.52.200
Nobel Brånemark®	CM - NB BM AURUM BASE (Desktop or Intraoral or on Ti base)		NP	41.14.001	41.50.001 41.52.001 41.53.100
					41.50.001 41.52.001 41.52.200
			RP	41.14.002	41.50.002 41.52.002 41.53.100
					41.50.002 41.52.002 41.52.200
			WP	41.14.003	41.50.003 41.52.003
					41.50.003 41.52.003
3i Osseotite®	CM - 3i OS AURUM BASE (Desktop or Intraoral or on Ti base)		NP/3.4	41.14.011	41.50.011 41.52.011 41.53.100
					41.50.011 41.52.011 41.52.200
			RP/4.1	41.14.012	41.50.012 41.52.012 41.53.100
					41.50.012 41.52.012 41.52.200
			WP/5.0	41.14.013	41.50.013 41.52.013 41.53.100
					41.50.013 41.52.013 41.52.200


 non-engaging


 non-engaging

 Scan Abutment
for desktop scanning
41.50.OXX


 Scan Abutment
for intra-oral scanning
41.52.OXX

















	Name Library	Subtype	Platform	Analogue	Scan Abutments		
3i Certain®	CM - 3i CE AURUM BASE (Desktop or Intraoral or on Ti base)		NP/3.4	41.14.014	41.50.014 41.52.014 41.53.100		
					41.50.014 41.52.014 41.52.200		
			RP/4.1	41.14.015	41.50.015 41.52.015 41.53.100		
					41.50.015 41.52.015 41.52.200		
			WP/5.0	41.14.016	41.50.016 41.52.016 41.53.100		
					41.50.016 41.52.016 41.52.200		
Straumann® Soft Tissue Level	CM - STR S AURUM BASE (Desktop or Intraoral or on Ti base)		RN/4.8	41.14.009	41.50.009 41.52.009		
					41.50.009 41.52.009		
			WN/6.5	41.14.045	41.50.045 41.52.045		
					41.50.045 41.52.045		
		Straumann® Bone Level	CM- STR BL AURUM BASE (Desktop or Intraoral or on Ti base)		NC/3.3	41.14.043	41.50.043 41.52.043 41.53.100
							41.50.043 41.52.043 41.52.200
	RC/4.1			41.14.044	41.50.044 41.52.044 41.53.100		
					41.50.044 41.52.044 41.52.200		
Zimmer Screw-Vent®	CM - ZI SCV AURUM BASE (Desktop or Intraoral or on Ti base)				NP/3.5	41.14.017	41.50.017 41.52.017 41.53.100
							41.50.017 41.52.017 41.52.200
			RP/4.5	41.14.018	41.50.018 41.52.018 41.53.100		
					41.50.018 41.52.018 41.52.200		
			WP/5.7	41.14.019	41.50.019 41.52.019 41.53.100		
					41.50.019 41.52.019 41.52.200		

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
41.50.0XX


 Scan Abutment
for intra-oral scanning
41.52.0XX


	Name Library	Type	Platform	Analogue	Scan Abutments
Astra tech implant system™ EV	CM - AS EV Ti-Base AURUM BASE (Desktop or Intraoral or on Ti base)		EV/3.6	41.14.059	41.50.059 41.52.059 41.53.100
					41.50.059 41.52.059 41.52.200
			EV/4.2	41.14.060	41.50.060 41.52.060 41.53.100
					41.50.060 41.52.060 41.52.200
			EV/4.8	41.14.061	41.50.061 41.52.061 41.53.100
					41.50.061 41.52.061 41.52.200
Astra tech osseospeed™	CM - AS AURUM BASE (Desktop or Intraoral or on Ti base)		RP	41.14.024	41.50.024 41.52.024 41.53.100
					41.50.024 41.52.024 41.52.200
			WP	41.14.025	41.50.025 41.52.025 41.53.100
					41.50.025 41.52.025 41.52.200
Friadent XIVE®	CM - DS XV AURUM BASE (Desktop or Intraoral or on Ti base)		NP	41.14.038	41.50.038 41.52.038 41.53.100
					41.50.038 41.52.038 41.52.200
			RP	41.14.039	41.50.039 41.52.039 41.53.100
					41.50.039 41.52.039 41.52.200
			WP	41.14.040	41.50.040 41.52.040 41.53.100
					41.50.040 41.52.040 41.52.200











 non-engaging


 non-engaging

 Scan Abutment
for desktop scanning
41.50.0XX

 Scan Abutment
for intra-oral scanning
41.52.0XX


 Scan Abutment
on Ti base
41.53.100 41.52.200

	Name Library	Type	Platform	Analogue	Scan Abutments
Camlog®	CM - CL AURUM BASE (Desktop or Intraoral or on Ti base)		3.8	14.065	50.065 52.065 53.100
					50.065 52.065 53.200
			4.3	14.066	50.066 52.066 53.100
					50.066 52.066 53.200
Mis® Seven	CM - MIS SVN AURUM BASE (Desktop or Intraoral or on Ti base)		SP	14.095	50.017 52.017 53.100
					50.017 52.017 53.200
			WP	14.096	50.018 52.018 53.100
					50.018 52.018 53.200
Multi-Unit	CM - Multi AURUM BASE (Desktop or Intraoral or on Ti base)		RP	14.007	50.007 52.007
			WP	14.047	50.047 52.047

 non-engaging

 non-engaging

 Scan Abutment
for desktop scanning
41.50.0XX

 Scan Abutment
for intra-oral scanning
41.52.0XX









 Scan Abutment
on Ti base
41.53.100 41.52.200

Table 1.5
CM Connect scan on Uniabutment &
Multi-Unit & Synocta

	Category	Type	Platform	Analogue	Scan Abutments
Uniabutment® on Astra Tech™ Osseospeed™	CM UNIABUTMENT DIRECT (Uniabutment Scan Abutment)		Uniabutment® 20°	41.14.036	41.50.036 41.52.036
	CM UNIABUTMENT DIRECT (Uniabutment Scan Abutment)		Uniabutment® 45°	41.14.037	41.50.037 41.52.037
Astra Tech™ Implant System™ EV	CM UNIABUTMENT DIRECT (Uniabutment Scan Abutment)		Uniabutment®	41.14.058	41.50.058 41.52.058
Multi-Unit	CM MULTIUNIT DIRECT (Multi-Unit Scan Abutment)		RP	41.14.007	41.50.036 41.52.036
	CM MULTIUNIT DIRECT (Multi-Unit Scan Abutment)		WP	41.14.047	41.50.037 41.52.037
Straumann® Soft Tissue level witch synOcta® Abutment	CM STR SYN DIRECT (Bars and bridges)		RN on synOcta®	41.14.010	41.50.010 41.52.010
	CM STR SYN DIRECT (Bars and bridges)		WN on synOcta®	41.14.046	41.50.046 41.52.037






Table 1.6
CM Connect Pre-milled Abutment Ti

	Name Library	Implant Type	Analogue	Scan Abutments	
Nobel Active®	CM - NB AC Pre-milled Ti (Desktop or Intraoral)	3.0	41.14.068	41.52.068	41.53.100
		NP	41.14.041	41.50.041	41.52.041
		RP	41.14.042	41.50.042	41.52.042
		WP	41.14.069	41.52.069	41.53.100
Nobel Replace Select™	CM - NB RP Pre-milled Ti (Desktop or Intraoral)	NP	41.14.004	41.50.004	41.52.004
		RP	41.14.005	41.50.005	41.52.005
		WP	41.14.006	41.50.006	41.52.006
		6.0	41.14.097	41.52.097	41.53.100
Nobel Brånemark®	CM - NB BM Pre-milled Ti (Desktop or Intraoral)	NP	41.14.001	41.50.001	41.52.001
		RP	41.14.002	41.50.002	41.52.002
		WP	41.14.003	41.50.003	41.52.003
3i Osseotite®	CM - 3i OS Pre-milled Ti (Desktop or Intraoral)	NP	41.14.011	41.50.011	41.52.011
		RP	41.14.012	41.50.012	41.52.012
		WP	41.14.013	41.50.013	41.52.013
3i Certain®	CM - 3i CE Pre-milled Ti (Desktop or Intraoral)	NP	41.14.014	41.50.014	41.52.014
		RP	41.14.015	41.50.015	41.52.015
		WP	41.14.016	41.50.016	41.52.016
Straumann® Soft Tissue Level	CM - STR S Pre-milled Ti (Desktop or Intraoral)	NNC	41.14.070	41.52.069	41.53.100
		RN	41.14.009	41.50.009	41.52.009
		WN	41.14.045	41.50.045	41.52.045




non-engaging




non-engaging



Scan Abutment for desktop scanning
41.50.0XX




Scan Abutment for intra-oral scanning
41.52.0XX



Scan Abutment on Ti base
41.53.100 41.52.200

	Name Library	Implant Type	Analogue	Scan Abutments	
Straumann® Bone Level	CM - STR BL Pre-milled Ti (Desktop or Intraoral)	NC	41.14.043	41.50.043	41.52.043 41.53.100
		RC	41.14.044	41.50.044	41.52.044 41.53.100
Zimmer Screw-Vent®	CM - ZI SCV Pre-milled Ti (Desktop or Intraoral)	NP	41.14.017	41.50.017	41.52.017 41.53.100
		RP	41.14.018	41.50.018	41.52.018 41.53.100
		WP	41.14.019	41.50.019	41.52.019 41.53.100
Astra tech implant system™ EV	CM - AS EV Pre-milled Ti (Desktop or Intraoral)	EV/3.6	41.14.059	41.50.059	41.52.059 41.53.100
		EV/4.2	41.14.060	41.50.060	41.52.060 41.53.100
		EV/4.8	41.14.061	41.50.061	41.52.061 41.53.100
Astra tech implant Osseospeed™	CM - AS OS Pre-milled Ti (Desktop or Intraoral)	RP	41.14.024	41.50.024	41.52.024 41.53.100
		WP	41.14.025	41.50.025	41.52.025 41.53.100
Friadent XIVE®	CM - XV Pre-milled Ti (Desktop or Intraoral)	NP	41.14.038	41.50.038	41.52.038 41.53.100
		RP	41.14.039	41.50.039	41.52.039 41.53.100
		WP	41.14.040	41.50.040	41.52.040 41.53.100
Ankylos® C/X	CM - ANK Pre-milled Ti (Desktop or Intraoral)	All platforms, height GH 2,5	41.14.048 41.14.049	41.50.049	41.52.049 41.53.100
Osstem® TS	CM - OS Pre-milled Ti (Desktop or Intraoral)	Mini Yellow	41.14.072	41.50.072	41.52.072 41.53.100
		STD Regular	41.14.073	41.50.073	41.52.073 41.53.100



non-engaging



non-engaging



Scan Abutment
for desktop scanning
41.50.0XX




Scan Abutment
for intra-oral scanning
41.52.0XX




Scan Abutment
on Ti base
41.53.100 41.52.200

	Name Library	Implant Type	Analogue	Scan Abutments
Camlog®	CM - CL Pre-milled Ti (Desktop or Intraoral)	3.3	41.14.064	41.52.064 41.53.100
		3.8	41.14.065	41.50.065 41.52.065 41.53.100
		4.3	41.14.066	41.50.066 41.52.066 41.53.100
		5.0	41.14.067	41.52.067 41.53.100
Biohorizns® External	CM - BH EXT Pre-milled Ti (Desktop or Intraoral)	RP/4.2	41.14.021	41.50.021 41.52.021 41.53.100
		WP/5.2	41.14.022	41.50.022 41.52.022
MIS® Seven	CM - MIS SVN Pre-milled Ti (Desktop or Intraoral)	SP	41.14.095	41.50.017 41.52.017 41.53.100
		WP	41.14.096	41.50.018 41.52.018 41.53.100
Megagen AntRidge®	CM - MG Pre-milled Ti (Desktop or Intraoral)		41.14.057	41.50.057 41.53.100
Dentium Superline™ & Implantium	CM - DT Pre-milled Ti (Desktop or Intraoral)	RP	41.14.090	41.52.090 41.53.100



non-engaging



non-engaging



Scan Abutment
for desktop scanning
41.50.OXX




Scan Abutment
for intra-oral scanning
41.52.OXX




Scan Abutment
on Ti base
41.53.100 41.52.200

Table 1.7
CM Connect Pre-milled Abutment CoCr

	Name Library	Implant Type	Analogue	Scan Abutments
Nobel Active®	CM - NB AC Pre-milled CoCr (Desktop or Intraoral)	3.0	41.14.068	41.52.068 41.53.100
		NP	41.14.041	41.50.041 41.52.041 41.53.100
		RP	41.14.042	41.50.042 41.52.042 41.53.100
		WP	41.14.069	41.52.069 41.53.100
Nobel Replace Select™	CM - NB RP Pre-milled CoCr (Desktop or Intraoral)	NP	41.14.004	41.50.004 41.52.004 41.53.100
		RP	41.14.005	41.50.005 41.52.005 41.53.100
		WP	41.14.006	41.50.006 41.52.006 41.53.100
		6.0	41.14.097	41.52.097 41.53.100
Nobel Brånemark®	CM - NB BM Pre-milled CoCr (Desktop or Intraoral)	NP	41.14.001	41.50.001 41.52.001 41.53.100
		RP	41.14.002	41.50.002 41.52.002 41.53.100
		WP	41.14.003	41.50.003 41.52.003
3i Osseotite®	CM - 3i OS Pre-milled CoCr (Desktop or Intraoral)	NP	41.14.011	41.50.011 41.52.011 41.53.100
		RP	41.14.012	41.50.012 41.52.012 41.53.100
		WP	41.14.013	41.50.013 41.52.013 41.53.100
3i Certain®	CM - 3i CE Pre-milled CoCr (Desktop or Intraoral)	NP	41.14.014	41.50.014 41.52.014 41.53.100
		RP	41.14.015	41.50.015 41.52.015 41.53.100
		WP	41.14.016	41.50.016 41.52.016 41.53.100
Straumann® Soft Tissue Level	CM - STR S Pre-milled CoCr (Desktop or Intraoral)	RN	41.14.009	41.50.009 41.52.009
		WN	41.14.045	41.50.045 41.52.045




non-engaging



non-engaging



Scan Abutment
for desktop scanning
41.50.0XX



Scan Abutment
for intra-oral scanning
41.52.0XX



Scan Abutment
on Ti base
41.53.100 41.52.200

	Name Library	Implant Type	Analogue	Scan Abutments	
Straumann® Bone Level	CM - STR BL Pre-milled CoCr (Desktop or Intraoral)	NC	41.14.043	41.50.043	41.52.043 41.53.100
		RC	41.14.044	41.50.044	41.52.044 41.53.100
Zimmer Screw-Vent®	CM - ZI SCV Pre-milled CoCr (Desktop or Intraoral)	NP	41.14.017	41.50.017	41.52.017 41.53.100
		RP	41.14.018	41.50.018	41.52.018 41.53.100
		WP	41.14.019	41.50.019	41.52.019 41.53.100
Astra tech implant system™ EV	CM - AS EV Pre-milled CoCr (Desktop or Intraoral)	EV/3.6	41.14.059	41.50.059	41.52.059 41.53.100
		EV/4.2	41.14.060	41.50.060	41.52.060 41.53.100
		EV/4.8	41.14.061	41.50.061	41.52.061 41.53.100
Astra tech implant Osseospeed™	CM - AS OS Pre-milled CoCr (Desktop or Intraoral)	RP	41.14.024	41.50.024	41.52.024 41.53.100
		WP	41.14.025	41.50.025	41.52.025 41.53.100
Friadent XIVE®	CM - XV Pre-milled CoCr (Desktop or Intraoral)	NP	41.14.038	41.50.038	41.52.038 41.53.100
		RP	41.14.039	41.50.039	41.52.039 41.53.100
		WP	41.14.040	41.50.040	41.52.040 41.53.100
Ankylos® C/X	CM - ANK Pre-milled CoCr (Desktop or Intraoral)	All platforms, height GH 2,5	41.14.048 41.14.049	41.50.049	41.52.049 41.53.100
		Mini Yellow	41.14.072	41.50.072	41.52.072 41.53.100
Osstem® TS	CM - OS Pre-milled CoCr (Desktop or Intraoral)	STD Regular	41.14.073	41.50.073	41.52.073 41.53.100



non-engaging



non-engaging



Scan Abutment
for desktop scanning
41.50.OXX




Scan Abutment
for intra-oral scanning
41.52.OXX



Scan Abutment
on Ti base
41.53.100 41.52.200

	Name Library	Implant Type	Analogue	Scan Abutments
Camlog®	CM - CL Pre-milled CoCr (Desktop or Intraoral)	3.3	41.14.064	41.52.064 41.53.100
		3.8	41.14.065	41.50.065 41.52.065 41.53.100
		4.3	41.14.066	41.50.066 41.52.066 41.53.100
		5.0	41.14.067	41.52.067 41.53.100
Biohorizns® External	CM - BH EXT Pre-milled CoCr (Desktop or Intraoral)	RP/4.2	41.14.021	41.50.021 41.52.021 41.53.100
		WP/5.2	41.14.022	41.50.022 41.52.022
MIS® Seven	CM - MIS SVN Pre-milled CoCr (Desktop or Intraoral)	SP	41.14.095	41.50.017 41.52.017 41.53.100
		WP	41.14.096	41.50.018 41.52.018 41.53.100
Megagen AntRidge®	CM - MG Pre-milled CoCr (Desktop or Intraoral)		41.14.057	41.50.057 41.53.100
Dentium Superline™ & Implantium	CM - DT Pre-milled CoCr (Desktop or Intraoral)	RP	41.14.090	41.52.090 41.53.100



non-engaging



non-engaging



Scan Abutment
for desktop scanning
41.50.OXX



Scan Abutment
for intra-oral scanning
41.52.OXX



Scan Abutment
on TI base
41.53.100 41.52.200

All CM Connect products are manufactured according to ISO 9001:2015 and ISO 13485:2016 with CE marking for all products Class IIb

All CM Connect AURUMBase® screws and screwdrivers are based on the standard ISO 10664

3i® and 3i Certain®, Miniplant®, MicroMiniplant™ and Osseotite® are registered trademarks of 3i Implant Innovations, USA
ASTRA TECH™, OsseoSpeed™, ANKYLOS® C/X, XiVE®, FRIADENT®, ASTRATECH implant system™ EV and UniAbutment®
are registered trademarks of DENTSPLY SIRONA Implants

BioHorizons® is a registered trademark of BioHorizons, INC.

Brånemark System®, NobelReplace® and Replace Select™, NobelActive™, NobelReplace® CC, NobelSpeedy™ and Multi-Unit are
registered trademarks of Nobel Biocare AB, Sweden

ITI®, Straumann® and synOcta® are registered trademarks of Straumann Holding AG, Switzerland

MIS® is a registered trademark of MIS Implant Technologies Ltd.

Zimmer Screw-Vent® and SwissPlus® are registered trademarks of Zimmer Dental INC, USA

Camlog® and Conelog® are registered trademarks of CAMLOG Biotechnologies AG

Megagen® is a registered trademark of MEGAGEN IMPLANT CO. LTD

TORX® is a trademark of ACUMENT INTELLECTUAL PROPERTIES, LLC

Bego Semados® is a registered trademark of BEGO GmbH & Co. KG

Anthogyr® is a registered trademark of Anthogyr SAS

Osstem® is a registered trademark of OSSTEM IMPLANT CO. LTD

Dentium® is a trademark of Dentium Co. Ltd.

Biotech Kontakt® is a trademark Biotech Dental SAS.

BTI® is a trademark BTI Biotechnology Institute.

CM Connect products can be used only by licensed dentists, physicians and labs.

Please refer ever to the instructions of use available on the website

www.cmsa.ch

or supplied with the product.

Some products may not be available in your country. Some of the products are cleared for US and Canada.

Please contact your local distributor for more details.