

## Description

Grinding points made of pure aluminium oxide for Palladium based and non precious bonding alloys.

## Indication

Abrasives for Palladium-base and non precious alloys with less than 50 % gold content. Surface grinding of crowns and bridges.

## Characteristics

- no negative effects due to carbon (diamond) remnants
- fine cutting grain structure
- wide range of cutting speed (5000–12000 rpm)
- true running due to precision centred shaft (∅ 2.35mm)
- long use and constant grain structure during wear

## Instructions for use

Work out the metal surface only in one direction with moderate pressure.

## Safety precautions

Do not breathe dust

Presentation  
box with 12 points  
08052193 Cerasiv blue

Rx only