



+ CM LOC®  
Implant compatibility list.

# CM LOC® Implantat Kompatibilitätsliste

Alle mit ® und ™ gekennzeichneten Unternehmen und Marken sind eingetragene Marken ihrer jeweiligen Inhaber.



## CM LOC® Systèmes d'implants compatibles

Toutes les sociétés et marques accompagnées des symboles ® et ™ sont des marques déposées de leurs propriétaires respectifs.

## CM LOC® Compatible implant systems

All companies and trademarks with ® and ™ are registered trademarks of their respective owners.

## CM LOC® Sistemi implantari compatibili

Tutte le ditte e i marchi contrassegnati con ® e ™ sono marchi registrati dei rispettivi proprietari.

## CM LOC® Sistemas de implantes compatibles

Todas las empresas y marcas comerciales designadas con ® y ™ son marcas registradas de sus correspondientes propietarios.

|            |              |            | Platform                    | Drehmoment<br>(empfohlen) /<br>Torque<br>(recommended) | Art. Nr.<br>Cendres + Métaux |
|------------|--------------|------------|-----------------------------|--|------------------------------|
| Straumann® | CM LOC®      | Abutment   | Straumann® RN 4.8 GH1       | 35Ncm  | 0500 1408                    |
|            | CM LOC®      | Abutment   | Straumann® RN 4.8 GH2       | 35Ncm  | 0500 1409                    |
|            | CM LOC®      | Abutment   | Straumann® RN 4.8 GH3       | 35Ncm  | 0500 1410                    |
|            | CM LOC®      | Abutment   | Straumann® RN 4.8 GH4       | 35Ncm  | 0500 1411                    |
|            | CM LOC®      | Abutment   | Straumann® RN 4.8 GH5       | 35Ncm  | 0500 1412                    |
|            | CM LOC®      | Case Guide | Straumann® RN 4.8           | –  | 0500 1578                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RN 4.8 GH1       | 35Ncm  | 0500 1720                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RN 4.8 GH2       | 35Ncm  | 0500 1721                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RN 4.8 GH3       | 35Ncm  | 0500 1722                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RN 4.8 GH4       | 35Ncm  | 0500 1723                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RN 4.8 GH5       | 35Ncm  | 0500 1724                    |
|            | CM LOC® FLEX | Case Guide | Straumann® RN 4.8           | –  | 0500 1884                    |
| Straumann® | CM LOC®      | Abutment   | Straumann® RC 4.1 / 4.8 GH1 | 35Ncm  | 0500 1413                    |
|            | CM LOC®      | Abutment   | Straumann® RC 4.1 / 4.8 GH2 | 35Ncm  | 0500 1414                    |
|            | CM LOC®      | Abutment   | Straumann® RC 4.1 / 4.8 GH3 | 35Ncm  | 0500 1415                    |
|            | CM LOC®      | Abutment   | Straumann® RC 4.1 / 4.8 GH4 | 35Ncm  | 0500 1416                    |
|            | CM LOC®      | Abutment   | Straumann® RC 4.1 / 4.8 GH5 | 35Ncm  | 0500 1417                    |
|            | CM LOC®      | Case Guide | Straumann® RC 4.1 / 4.8     | –  | 0500 1579                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RC 4.1 / 4.8 GH1 | 35Ncm  | 0500 1725                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RC 4.1 / 4.8 GH2 | 35Ncm  | 0500 1726                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RC 4.1 / 4.8 GH3 | 35Ncm  | 0500 1727                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RC 4.1 / 4.8 GH4 | 35Ncm  | 0500 1728                    |
|            | CM LOC® FLEX | Abutment   | Straumann® RC 4.1 / 4.8 GH5 | 35Ncm  | 0500 1729                    |
|            | CM LOC® FLEX | Case Guide | Straumann® RC 4.1 / 4.8     | –  | 0500 1885                    |
| Straumann® | CM LOC®      | Abutment   | Straumann® NNC 3.5 GH1      | 35Ncm  | 0500 1418                    |
|            | CM LOC®      | Abutment   | Straumann® NNC 3.5 GH2      | 35Ncm  | 0500 1419                    |
|            | CM LOC®      | Abutment   | Straumann® NNC 3.5 GH3      | 35Ncm  | 0500 1420                    |
|            | CM LOC®      | Abutment   | Straumann® NNC 3.5 GH4      | 35Ncm  | 0500 1421                    |
|            | CM LOC®      | Abutment   | Straumann® NNC 3.5 GH5      | 35Ncm  | 0500 1422                    |
|            | CM LOC®      | Case Guide | Straumann® NNC 3.5          | –  | 0500 1580                    |
| Straumann® | CM LOC®      | Abutment   | Straumann® NC 3.3 GH1       | 35Ncm  | 0500 2547                    |
|            | CM LOC®      | Abutment   | Straumann® NC 3.3 GH2       | 35Ncm  | 0500 2548                    |
|            | CM LOC®      | Abutment   | Straumann® NC 3.3 GH3       | 35Ncm  | 0500 2549                    |
|            | CM LOC®      | Abutment   | Straumann® NC 3.3 GH4       | 35Ncm  | 0500 2550                    |
|            | CM LOC®      | Abutment   | Straumann® NC 3.3 GH5       | 35Ncm  | 0500 2551                    |
|            | CM LOC®      | Case Guide | Straumann® NC 3.3           | –  | 0500 2552                    |

|               |               |            | Plattform                                | Drehmoment<br>(empfohlen) /<br>Torque<br>(recommended) | Art. Nr.<br>Cendres + Métaux |
|---------------|---------------|------------|--|--|------------------------------|
| Nobel Biocare | CM LOC®       | Abutment   | Nobel Biocare Replace Select® NP 3.5 GH1 | 35Ncm  | 0500 1423                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Replace Select® NP 3.5 GH2 | 35Ncm  | 0500 1424                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Replace Select® NP 3.5 GH3 | 35Ncm  | 0500 1425                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Replace Select® NP 3.5 GH4 | 35Ncm  | 0500 1426                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Replace Select® NP 3.5 GH5 | 35Ncm  | 0500 1427                    |
|               | CM LOC®       | Case Guide | Nobel Biocare Replace Select® NP 3.5     | –  | 0500 1581                    |
| Nobel Biocare | CM LOC®       | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH1 | 35Ncm  | 0500 1296                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH2 | 35Ncm  | 0500 1300                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH3 | 35Ncm  | 0500 1301                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH4 | 35Ncm  | 0500 1302                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH5 | 35Ncm  | 0500 1303                    |
|               | CM LOC®       | Case Guide | Nobel Biocare Replace Select® RP 4.3     | –  | 0500 1582                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH1 | 35Ncm  | 0500 1740                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH2 | 35Ncm  | 0500 1741                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH3 | 35Ncm  | 0500 1742                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH4 | 35Ncm  | 0500 1743                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Replace Select® RP 4.3 GH5 | 35Ncm  | 0500 1744                    |
|               | CM LOC® FLEX  | Case Guide | Nobel Biocare Replace Select® RP 4.3     | –  | 0500 1888                    |
|               | Nobel Biocare | CM LOC®    | Abutment                                 | Nobel Biocare Active® NP 3.5 GH1                       | 35Ncm                        |
| CM LOC®       |               | Abutment   | Nobel Biocare Active® NP 3.5 GH2         | 35Ncm  | 0500 2110                    |
| CM LOC®       |               | Abutment   | Nobel Biocare Active® NP 3.5 GH3         | 35Ncm  | 0500 2111                    |
| CM LOC®       |               | Abutment   | Nobel Biocare Active® NP 3.5 GH4         | 35Ncm  | 0500 2112                    |
| CM LOC®       |               | Abutment   | Nobel Biocare Active® NP 3.5 GH5         | 35Ncm  | 0500 2113                    |
| CM LOC®       |               | Case Guide | Nobel Biocare Active® NP 3.5             | –  | 0500 2187                    |
| Nobel Biocare | CM LOC®       | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH1     | 35Ncm  | 0500 1437                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH2     | 35Ncm  | 0500 1438                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH3     | 35Ncm  | 0500 1439                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH4     | 35Ncm  | 0500 1440                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH5     | 35Ncm  | 0500 1441                    |
|               | CM LOC®       | Case Guide | Nobel Biocare Active® RP 4.3/5.0         | –  | 0500 1584                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH1     | 35Ncm  | 0500 1750                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH2     | 35Ncm  | 0500 1751                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH3     | 35Ncm  | 0500 1752                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH4     | 35Ncm  | 0500 1753                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Active® RP 4.3/5.0 GH5     | 35Ncm  | 0500 1754                    |
|               | CM LOC® FLEX  | Case Guide | Nobel Biocare Active® RP 4.3/5.0         | –  | 0500 1890                    |
| Nobel Biocare | CM LOC®       | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH1      | 35Ncm  | 0500 2114                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH2      | 35Ncm  | 0500 2115                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH3      | 35Ncm  | 0500 2116                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH4      | 35Ncm  | 0500 2117                    |
|               | CM LOC®       | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH5      | 35Ncm  | 0500 2118                    |
|               | CM LOC®       | Case Guide | Nobel Biocare Brånemark® RP 4.0          | –  | 0500 2188                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH1      | 35Ncm  | 0500 2261                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH2      | 35Ncm  | 0500 2262                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH3      | 35Ncm  | 0500 2263                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH4      | 35Ncm  | 0500 2264                    |
|               | CM LOC® FLEX  | Abutment   | Nobel Biocare Brånemark® RP 4.0 GH5      | 35Ncm  | 0500 2265                    |
|               | CM LOC® FLEX  | Case Guide | Nobel Biocare Brånemark® RP 4.0          | –  | 0500 2309                    |

|            |              |            | Platform                           | Drehmoment<br>(empfohlen) /<br>Torque<br>(recommended) | Art. Nr.<br>Cendres + Métaux |
|------------|--------------|------------|------------------------------------|--|------------------------------|
| Astra Tech | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 3.5/4.0 GH1 | 25Ncm  | 0500 1452                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 3.5/4.0 GH2 | 25Ncm  | 0500 1453                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 3.5/4.0 GH3 | 25Ncm  | 0500 1454                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 3.5/4.0 GH4 | 25Ncm  | 0500 1455                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 3.5/4.0 GH5 | 25Ncm  | 0500 1456                    |
|            | CM LOC®      | Case Guide | Astra Tech OsseoSpeed® 3.5/4.0     | –  | 0500 1587                    |
| Astra Tech | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH1 | 25Ncm  | 0500 1457                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH2 | 25Ncm  | 0500 1458                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH3 | 25Ncm  | 0500 1459                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH4 | 25Ncm  | 0500 1460                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH5 | 25Ncm  | 0500 1461                    |
|            | CM LOC®      | Case Guide | Astra Tech OsseoSpeed® 4.5/5.0     | –  | 0500 1588                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH1 | 25Ncm  | 0500 1770                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH2 | 25Ncm  | 0500 1771                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH3 | 25Ncm  | 0500 1772                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH4 | 25Ncm  | 0500 1773                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® 4.5/5.0 GH5 | 25Ncm  | 0500 1774                    |
|            | CM LOC® FLEX | Case Guide | Astra Tech OsseoSpeed® 4.5/5.0     | –  | 0500 1894                    |
| Astra Tech | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 3.6 GH1  | 25Ncm  | 0500 1963                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 3.6 GH2  | 25Ncm  | 0500 1964                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 3.6 GH3  | 25Ncm  | 0500 1965                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 3.6 GH4  | 25Ncm  | 0500 1966                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 3.6 GH5  | 25Ncm  | 0500 1967                    |
|            | CM LOC®      | Case Guide | Astra Tech OsseoSpeed® EV 3.6      | –  | 0500 1981                    |
| Astra Tech | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH1  | 25Ncm  | 0500 1968                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH2  | 25Ncm  | 0500 1969                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH3  | 25Ncm  | 0500 1970                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH4  | 25Ncm  | 0500 1971                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH5  | 25Ncm  | 0500 1972                    |
|            | CM LOC®      | Case Guide | Astra Tech OsseoSpeed® EV 4.2      | –  | 0500 1982                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH1  | 25Ncm  | 0500 1932                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH2  | 25Ncm  | 0500 1933                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH3  | 25Ncm  | 0500 1934                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH4  | 25Ncm  | 0500 1935                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.2 GH5  | 25Ncm  | 0500 1936                    |
|            | CM LOC® FLEX | Case Guide | Astra Tech OsseoSpeed® EV 4.2      | –  | 0500 1946                    |
| Astra Tech | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH1  | 25Ncm  | 0500 1973                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH2  | 25Ncm  | 0500 1974                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH3  | 25Ncm  | 0500 1975                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH4  | 25Ncm  | 0500 1976                    |
|            | CM LOC®      | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH5  | 25Ncm  | 0500 1977                    |
|            | CM LOC®      | Case Guide | Astra Tech OsseoSpeed® EV 4.8      | –  | 0500 1983                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH1  | 25Ncm  | 0500 1937                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH2  | 25Ncm  | 0500 1938                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH3  | 25Ncm  | 0500 1939                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH4  | 25Ncm  | 0500 1940                    |
|            | CM LOC® FLEX | Abutment   | Astra Tech OsseoSpeed® EV 4.8 GH5  | 25Ncm  | 0500 1941                    |
|            | CM LOC® FLEX | Case Guide | Astra Tech OsseoSpeed® EV 4.8      | –  | 0500 1947                    |

|   |              |            | Plattform  | Drehmoment<br>(empfohlen) /<br>Torque<br>(recommended) | Art. Nr.<br>Cendres + Métaux |
|---|--------------|------------|--|--|------------------------------|
| Zimmer<br>MIS® (standard<br>narrow)<br>BioHorizons® | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 3.5 GH1 | 30Ncm  | 0500 1462                    |
|   | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 3.5 GH2 | 30Ncm  | 0500 1463                    |
|   | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 3.5 GH3 | 30Ncm  | 0500 1464                    |
|   | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 3.5 GH4 | 30Ncm  | 0500 1465                    |
|   | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 3.5 GH5 | 30Ncm  | 0500 1466                    |
|   | CM LOC®      | Case Guide | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 3.5     | –  | 0500 1589                    |
| Zimmer<br>MIS® (wide plat-<br>form)<br>BioHorizons® | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH1 | 30Ncm  | 0500 1467                    |
|   | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH2 | 30Ncm  | 0500 1468                    |
|   | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH3 | 30Ncm  | 0500 1469                    |
|   | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH4 | 30Ncm  | 0500 1470                    |
|   | CM LOC®      | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH5 | 30Ncm  | 0500 1471                    |
|   | CM LOC®      | Case Guide | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5     | –  | 0500 1590                    |
|   | CM LOC® FLEX | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH1 | 30Ncm  | 0500 1780                    |
|   | CM LOC® FLEX | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH2 | 30Ncm  | 0500 1781                    |
|   | CM LOC® FLEX | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH3 | 30Ncm  | 0500 1782                    |
|   | CM LOC® FLEX | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH4 | 30Ncm  | 0500 1783                    |
|   | CM LOC® FLEX | Abutment   | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5 GH5 | 30Ncm  | 0500 1784                    |
|   | CM LOC® FLEX | Case Guide | Zimmer Tapered Screw-Vent® / MiS® Seven / BioHorizons® 4.5     | –  | 0500 1896                    |

|          |              |            | Platform                                   | Drehmoment<br>(empfohlen) /<br>Torque<br>(recommended) | Art. Nr.<br>Cendres + Métaux |
|----------|--------------|------------|--|--|------------------------------|
| Camlog®  | CM LOC®      | Abutment   | Camlog® 3.8 GH1                            | 30Ncm  | 0500 1544                    |
|          | CM LOC®      | Abutment   | Camlog® 3.8 GH2                            | 30Ncm  | 0500 1545                    |
|          | CM LOC®      | Abutment   | Camlog® 3.8 GH3                            | 30Ncm  | 0500 1546                    |
|          | CM LOC®      | Abutment   | Camlog® 3.8 GH4                            | 30Ncm  | 0500 1547                    |
|          | CM LOC®      | Case Guide | Camlog® 3.8                                | –  | 0500 1591                    |
|          | CM LOC® FLEX | Abutment   | Camlog® 3.8 GH1                            | 30Ncm  | 0500 1790                    |
|          | CM LOC® FLEX | Abutment   | Camlog® 3.8 GH2                            | 30Ncm  | 0500 1791                    |
|          | CM LOC® FLEX | Abutment   | Camlog® 3.8 GH3                            | 30Ncm  | 0500 1792                    |
|          | CM LOC® FLEX | Abutment   | Camlog® 3.8 GH4                            | 30Ncm  | 0500 1793                    |
|          | CM LOC® FLEX | Case Guide | Camlog® 3.8                                | –  | 0500 1897                    |
| Camlog®  | CM LOC®      | Abutment   | Camlog® 4.3 GH1                            | 30Ncm  | 0500 1549                    |
|          | CM LOC®      | Abutment   | Camlog® 4.3 GH2                            | 30Ncm  | 0500 1550                    |
|          | CM LOC®      | Abutment   | Camlog® 4.3 GH3                            | 30Ncm  | 0500 1551                    |
|          | CM LOC®      | Abutment   | Camlog® 4.3 GH4                            | 30Ncm  | 0500 1552                    |
|          | CM LOC®      | Case Guide | Camlog® 4.3                                | –  | 0500 1592                    |
|          | CM LOC® FLEX | Abutment   | Camlog® 4.3 GH1                            | 30Ncm  | 0500 1795                    |
|          | CM LOC® FLEX | Abutment   | Camlog® 4.3 GH2                            | 30Ncm  | 0500 1796                    |
|          | CM LOC® FLEX | Abutment   | Camlog® 4.3 GH3                            | 30Ncm  | 0500 1797                    |
|          | CM LOC® FLEX | Abutment   | Camlog® 4.3 GH4                            | 30Ncm  | 0500 1798                    |
|          | CM LOC® FLEX | Case Guide | Camlog® 4.3                                | –  | 0500 1898                    |
| Camlog®  | CM LOC®      | Abutment   | Conelog® 3.8/4.3 GH1                       | 30Ncm  | 0500 1482                    |
|          | CM LOC®      | Abutment   | Conelog® 3.8/4.3 GH2                       | 30Ncm  | 0500 1483                    |
|          | CM LOC®      | Abutment   | Conelog® 3.8/4.3 GH3                       | 30Ncm  | 0500 1484                    |
|          | CM LOC®      | Abutment   | Conelog® 3.8/4.3 GH4                       | 30Ncm  | 0500 1485                    |
|          | CM LOC®      | Abutment   | Conelog® 3.8/4.3 GH5                       | 30Ncm  | 0500 1486                    |
|          | CM LOC®      | Case Guide | Conelog® 3.8/4.3                           | –  | 0500 1594                    |
|          | CM LOC® FLEX | Abutment   | Conelog® 3.8/4.3 GH1                       | 30Ncm  | 0500 1805                    |
|          | CM LOC® FLEX | Abutment   | Conelog® 3.8/4.3 GH2                       | 30Ncm  | 0500 1806                    |
|          | CM LOC® FLEX | Abutment   | Conelog® 3.8/4.3 GH3                       | 30Ncm  | 0500 1807                    |
|          | CM LOC® FLEX | Abutment   | Conelog® 3.8/4.3 GH4                       | 30Ncm  | 0500 1808                    |
|          | CM LOC® FLEX | Abutment   | Conelog® 3.8/4.3 GH5                       | 30Ncm  | 0500 1809                    |
|          | CM LOC® FLEX | Case Guide | Conelog® 3.8/4.3                           | –  | 0500 1899                    |
| Dentsply | CM LOC®      | Abutment   | Dentsply Ankylos® C 3.5, 4.5, 5.5, 7.0 GH1 | 25Ncm  | 0500 2014                    |
|          | CM LOC®      | Abutment   | Dentsply Ankylos® C 3.5, 4.5, 5.5, 7.0 GH2 | 25Ncm  | 0500 2015                    |
|          | CM LOC®      | Abutment   | Dentsply Ankylos® C 3.5, 4.5, 5.5, 7.0 GH3 | 25Ncm  | 0500 2016                    |
|          | CM LOC®      | Abutment   | Dentsply Ankylos® C 3.5, 4.5, 5.5, 7.0 GH4 | 25Ncm  | 0500 2017                    |
|          | CM LOC®      | Abutment   | Dentsply Ankylos® C 3.5, 4.5, 5.5, 7.0 GH5 | 25Ncm  | 0500 2018                    |
|          | CM LOC®      | Case Guide | Dentsply Ankylos® C 3.5, 4.5, 5.5, 7.0     | –  | 0500 2020                    |
|          | CM LOC® FLEX | Abutment   | Dentsply Ankylos® C 4.5, 5.5, 7.0 GH1      | 25Ncm  | 0500 2002                    |
|          | CM LOC® FLEX | Abutment   | Dentsply Ankylos® C 4.5, 5.5, 7.0 GH2      | 25Ncm  | 0500 2003                    |
|          | CM LOC® FLEX | Abutment   | Dentsply Ankylos® C 4.5, 5.5, 7.0 GH3      | 25Ncm  | 0500 2004                    |
|          | CM LOC® FLEX | Abutment   | Dentsply Ankylos® C 4.5, 5.5, 7.0 GH4      | 25Ncm  | 0500 2005                    |
|          | CM LOC® FLEX | Abutment   | Dentsply Ankylos® C 4.5, 5.5, 7.0 GH5      | 25Ncm  | 0500 2006                    |
|          | CM LOC® FLEX | Case Guide | Dentsply Ankylos® C 4.5, 5.5, 7.0          | –  | 0500 2008                    |

|                 |                 |            | Plattform                                     | Drehmoment<br>(empfohlen) /<br>Torque<br>(recommended) | Art. Nr.<br>Cendres + Métaux |
|-----------------|-----------------|------------|---|--|------------------------------|
| Sweden+ Martina | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.3 GH1         | 30Ncm  | 05002119                     |
|                 | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.3 GH2         | 30Ncm  | 05002120                     |
|                 | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.3 GH3         | 30Ncm  | 05002121                     |
|                 | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.3 GH4         | 30Ncm  | 05002122                     |
|                 | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.3 GH5         | 30Ncm  | 05002123                     |
|                 | CM LOC®         | Case Guide | Sweden+ Martina Premium Kohno 3.3             | –  | 05002189                     |
| Sweden+ Martina | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH1         | 30Ncm  | 05002124                     |
|                 | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH2         | 30Ncm  | 05002125                     |
|                 | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH3         | 30Ncm  | 05002126                     |
|                 | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH4         | 30Ncm  | 05002127                     |
|                 | CM LOC®         | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH5         | 30Ncm  | 05002128                     |
|                 | CM LOC®         | Case Guide | Sweden+ Martina Premium Kohno 3.8             | –  | 05002190                     |
|                 | CM LOC® FLEX    | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH1         | 30Ncm  | 05002266                     |
|                 | CM LOC® FLEX    | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH2         | 30Ncm  | 05002267                     |
|                 | CM LOC® FLEX    | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH3         | 30Ncm  | 05002268                     |
|                 | CM LOC® FLEX    | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH4         | 30Ncm  | 05002269                     |
|                 | CM LOC® FLEX    | Abutment   | Sweden+ Martina Premium Kohno 3.8 GH5         | 30Ncm  | 05002270                     |
|                 | CM LOC® FLEX    | Case Guide | Sweden+ Martina Premium Kohno 3.8             | –  | 05002310                     |
|                 | Sweden+ Martina | CM LOC®    | Abutment                                      | Sweden+ Martina Premium Kohno<br>4.25 GH1              | 30Ncm                        |
| CM LOC®         |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH2     | 30Ncm  | 05002130                     |
| CM LOC®         |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH3     | 30Ncm  | 05002131                     |
| CM LOC®         |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH4     | 30Ncm  | 05002132                     |
| CM LOC®         |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH5     | 30Ncm  | 05002133                     |
| CM LOC®         |                 | Case Guide | Sweden+ Martina Premium Kohno<br>4.25/5.0/6.0 | –  | 05002192                     |
| CM LOC® FLEX    |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH1     | 30Ncm  | 05002271                     |
| CM LOC® FLEX    |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH2     | 30Ncm  | 05002272                     |
| CM LOC® FLEX    |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH3     | 30Ncm  | 05002273                     |
| CM LOC® FLEX    |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH4     | 30Ncm  | 05002274                     |
| CM LOC® FLEX    |                 | Abutment   | Sweden+ Martina Premium Kohno<br>4.25 GH5     | 30Ncm  | 05002275                     |

|                  |              |            | Platform                                    | Drehmoment<br>(empfohlen) /<br>Torque<br>(recommended) | Art. Nr.<br>Cendres + Métaux |
|------------------|--------------|------------|---|--|------------------------------|
| Sweden + Martina | CM LOC®      | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH1  | 30Ncm  | 0500 2134                    |
|                  | CM LOC®      | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH2  | 30Ncm  | 0500 2135                    |
|                  | CM LOC®      | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH3  | 30Ncm  | 0500 2136                    |
|                  | CM LOC®      | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH4  | 30Ncm  | 0500 2137                    |
|                  | CM LOC®      | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH5  | 30Ncm  | 0500 2138                    |
|                  | CM LOC®      | Case Guide | Sweden + Martina Premium Kohno 5.0/6.0      | –  | 0500 2192                    |
|                  | CM LOC® FLEX | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH1  | 30Ncm  | 0500 2276                    |
|                  | CM LOC® FLEX | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH2  | 30Ncm  | 0500 2277                    |
|                  | CM LOC® FLEX | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH3  | 30Ncm  | 0500 2278                    |
|                  | CM LOC® FLEX | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH4  | 30Ncm  | 0500 2279                    |
|                  | CM LOC® FLEX | Abutment   | Sweden + Martina Premium Kohno 5.0/6.0 GH5  | 30Ncm  | 0500 2280                    |
|                  | CM LOC® FLEX | Case Guide | Sweden + Martina Premium Kohno 4.25/5.0/6.0 | –  | 0500 2312                    |
|                  | Osstem®      | CM LOC®    | Abutment                                    | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH1             | 30Ncm                        |
| CM LOC®          |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH2  | 30Ncm  | 0500 2160                    |
| CM LOC®          |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH3  | 30Ncm  | 0500 2161                    |
| CM LOC®          |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH4  | 30Ncm  | 0500 2162                    |
| CM LOC®          |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH5  | 30Ncm  | 0500 2163                    |
| CM LOC®          |              | Case Guide | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0      | –  | 0500 2197                    |
| CM LOC® FLEX     |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH1  | 30Ncm  | 0500 2566                    |
| CM LOC® FLEX     |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH2  | 30Ncm  | 0500 2567                    |
| CM LOC® FLEX     |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH3  | 30Ncm  | 0500 2568                    |
| CM LOC® FLEX     |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH4  | 30Ncm  | 0500 2569                    |
| CM LOC® FLEX     |              | Abutment   | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0 GH5  | 30Ncm  | 0500 2570                    |
| CM LOC® FLEX     |              | Case Guide | Osstem® TS Regular 4.0/4.5/5.0/6.0/7.0      | –  | 0500 2586                    |
| Osstem®          |              | CM LOC®    | Abutment                                    | Osstem® TS Mini 3.5 GH1                                | 30Ncm                        |
|                  | CM LOC®      | Abutment   | Osstem® TS Mini 3.5 GH2                     | 30Ncm  | 0500 2591                    |
|                  | CM LOC®      | Abutment   | Osstem® TS Mini 3.5 GH3                     | 30Ncm  | 0500 2592                    |
|                  | CM LOC®      | Abutment   | Osstem® TS Mini 3.5 GH4                     | 30Ncm  | 0500 2593                    |
|                  | CM LOC®      | Abutment   | Osstem® TS Mini 3.5 GH5                     | 30Ncm  | 0500 2594                    |
|                  | CM LOC®      | Case Guide | Osstem® TS Mini 3.5                         | –  | 0500 2596                    |
| Osstem®          | CM LOC®      | Abutment   | Osstem® US Regular 4.1 GH1                  | 30Ncm  | 0500 2154                    |
|                  | CM LOC®      | Abutment   | Osstem® US Regular 4.1 GH2                  | 30Ncm  | 0500 2155                    |
|                  | CM LOC®      | Abutment   | Osstem® US Regular 4.1 GH3                  | 30Ncm  | 0500 2156                    |
|                  | CM LOC®      | Abutment   | Osstem® US Regular 4.1 GH4                  | 30Ncm  | 0500 2157                    |
|                  | CM LOC®      | Abutment   | Osstem® US Regular 4.1 GH5                  | 30Ncm  | 0500 2158                    |
|                  | CM LOC®      | Case Guide | Osstem® US Regular 4.1                      | –  | 0500 2196                    |
| Osstem®          | CM LOC®      | Abutment   | Osstem® SS Regular 4.8 GH1                  | 30Ncm  | 0500 2149                    |
|                  | CM LOC®      | Abutment   | Osstem® SS Regular 4.8 GH2                  | 30Ncm  | 0500 2150                    |
|                  | CM LOC®      | Abutment   | Osstem® SS Regular 4.8 GH3                  | 30Ncm  | 0500 2151                    |
|                  | CM LOC®      | Abutment   | Osstem® SS Regular 4.8 GH4                  | 30Ncm  | 0500 2152                    |
|                  | CM LOC®      | Abutment   | Osstem® SS Regular 4.8 GH5                  | 30Ncm  | 0500 2153                    |
|                  | CM LOC®      | Case Guide | Osstem® SS Regular 4.8                      | –  | 0500 2195                    |



|         |                         | Platform  | Drehmoment<br>(empfohlen) /<br>Torque<br>(recommended) | Art. Nr.<br>Cendres + Métaux |
|---------|-------------------------|---|--|------------------------------|
| Dyna    | CM LOC®                 | for Dyna Helix Implant System<br>Distributed by Dyna Dental Engineering b.v,<br>The Netherlands |  |                              |
| Aurosan | CM LOC®<br>CM LOC® FLEX | for Aurosan REVOIS® PRO Implant System<br>Distributed by Aurosan GmbH, Germany                  |  |                              |