



**Precious metal milling.**  
Confidence meets technology –  
the best of two worlds.

# Precious metal milling.

The fusion of two core competences. CM expands the CAD/CAM milling service and extends the material portfolio.

## Confidence meets technology – the best of two worlds.



- easy realisation of complex superstructures in precious metal
- available on all major implants systems
- evidence based material for reconstructive prosthetic solutions




- high precision and stress-free fit thanks to CAD/CAM technology
- complex structures easy made

## Workflow.



By mail: [dental-cadcam@cmsa.ch](mailto:dental-cadcam@cmsa.ch)  
Data up: [www.cmsa.ch/cadcam-order](http://www.cmsa.ch/cadcam-order)  
By post: Cendres+ Métaux SA  
Scan+Design Centre  
Rue de Boujean 122  
P.O. Box  
CH-2501 Biel/Bienne

## Economical.



Thanks to the CAD/CAM technology the laboratory and the patient can save money.

- Significant material savings possible. Approximately 30% less material compared to the traditional process.
- Significant savings in laboratory working time, no more modelling, investing and casting procedures.
- Less use of consumables
- Easy processing, no more error prone procedures

## Indications.



- On Implants: crown, bridge, custom abutment
- On natural teeth: crown, bridge

## Compatibility.



Implant supported solutions are available on all major implant systems. Compatibility list according to:  
[www.cmsa.ch/en/Dental/Download-Center](http://www.cmsa.ch/en/Dental/Download-Center)

Alloy	Type	Au	Pt	Pd	Ag	Cu	Sn	Zn	In	Ga	Ir	Ru	Rh	Fe	CET (25–500°C) 10 <sup>-6</sup> K <sup>-1</sup>
Estetikor® Lumina PF	4	84.50	13.30					1.90			0.10		0.10	0.10	14.2
Estetikor® CC	4	12.00		52.29	23.00		2.00		10.00	0.50		0.20			14.2
Dentalor® 60	4	60.00	0.45	3.00	22.50	12.50		1.50			0.05				

## Contact.

Thanks to our local companies and wide distribution network, we are close by and available at all times.

[www.cmsa.ch/en/Dental/Contact](http://www.cmsa.ch/en/Dental/Contact)  
[www.cmsa.ch/en/Dental/Contact/Contact-form](http://www.cmsa.ch/en/Dental/Contact/Contact-form)  
[www.cmsa.ch/en/Dental/Contact/Order-form](http://www.cmsa.ch/en/Dental/Contact/Order-form)

Ask your local sales representative, he will be of any assistance.